

Technical Specifications

Model	OP830	OP330
Table Length (mm)	2040	2040
Table Width (w/o side rails, mm)	520	520
Table Height (w/o pads, mm)	680 - 950	680 - 950
Trendelenburg/Reverse	30°/30°	30°/30°
Lateral Tilt(left/right)	20°/20°	20°/20°
BackPlate Up/Down	75°/30°	75°/30°
Head Plate Up/Down (detachable)	30°/45°	30°/45°
Leg Plate Up/Down (detachable)	5°/90°	5°/90°
Leg Plate Split(detachable)	180°	180°
Flex/Reflex	210° / 105°	210° / 105°
Kidney-bridge Elevator Up (mm)	136	136
Longitudinal Shift (mm)	300	N/A
Weight Capacity (kg)	297	297



Corporate Headquarter
 HEYER Medical AG
 Carl-Heyer-Str. 1/3, 56130 Bad Ems, Germany
 Tel: +49 2603 791 0
 Fax: +49 2603 791 200
 info@heyemed.com
 www.heyemedical.de

Manufacturer
 HEYER Medical Co., Ltd.
 No.10, Chaobai Street, Yanjiao Development Zone
 065201, Sanhe City, Hebei, China
 Tel: +86 10 5841 1198
 Fax: +86 10 6371 8989
 marketing@heyemed.com

AOP830/330-2024.10



Aegistab OP830/330 Operating Table



A cost-effective electronic-hydraulic operating table with a safe load of 297kg, meeting the requirements of different surgeries.



① Premium Quality Material

- Table frame is specially coated to increase corrosion resistance
- 304 stainless steel is used for exposed core components, such as the central column and table base.



② Perfect C-arm Access

- 300mm longitudinal shifting capability creates perfect C-arm access




③ Ergonomic Adjustment

- Convenient leg plate and head plate positioning ensures one-hand control operation
- Imported gas springs guarantee smooth movement and accurate positioning
- Wide range of table height adjustment and T-shape base provide comfortable working environment



④ Surgical Configuration

- One-click flex and reflex positions
- Flexible patient positionings with optional accessories are suitable for a wide variety of surgeries



⑤ Intelligent safety system

- Remote controller panel will power off automatically when there is no operation over 1 minute
- Misoperation preventing design guarantees patients' safety during surgeries



⑥ Pad

- PU foam mattress, easy to clean and disinfection
- High-density material effectively reduces the possibilities of patient bedsores



⑦ Built-in Elevator Kidney-bridge

- Built-in kidney-bridge supports kidney and gall bladder surgeries efficiently (Optional X-ray radiolucent)