








Safe, Simple, Smart hawk-s1 Syringe pump




Simple

-  Intuitive interface with colorful screen, all-in-one information at a glance.
-  Scientific button design facilitates fast programming
Compact and light weight design allow easy moving and transferring.
-  Detachable handle and rotatable pole clamp.

Versatile

-  Automatically identifies various syringe types and sizes from 5ml to 60ml .
-  Pumps stack and lock together and can be expanded from 2 to 12.
-  Pressure and infusion process can be detected synchronously .
-  Illumination mode facilitates medical care at night.

Safe

-  Alarm management: With audible and visual alarm, caregivers can easily locate the alarm and warning.
-  Deliver the precise medication dosage to even a neonate or an infant.
Flow rates deliver with +/- 2% accuracy as low as 0.1 mL/hr.
-  Drug library: more than 1000 drugs in storage and recommended dose limits.

Simplified workflow

Wireless technology allows for connection to infusion management system HK-M1000 which features complete interoperability with HIS and CIS.



As Simple as it is Precise HK-400 Series Syringe Pump

- Intuitive interface with colorful screen, all-in-one information at a glance.
- Numeric keys combines button design facilitates fast programming.
- Compact and light weight design allow easy moving and transferring.
- Caregivers can easily install or remove pumps thanks to the integrated handle and rotatable pole clamp.
- Automatically identifies various syringe types and sizes from 5ml to 60ml .
- Numeric and graphic pressure indicator helps caregivers identify occlusion in advance to prevent potential risks.
- Alarm management: With audible and visual alarm, caregivers can easily locate the alarm and warning.
- Deliver the precise medication dosage to even your tiniest patients. Flow rates deliver with +/- 2% accuracy as low as 0.1 mL/hr.
- Long-lasting battery: the battery life is long-lasting enough to last up to 6 hours.
- Advanced moisture protection(IP24)for cleaning and disinfection results in optimized stability and fewer repair costs.
- Wireless technology allows for connection to infusion management system HK-M1000 which features complete interoperability with HIS and CIS.

